



Formalin is among other things **carcinogenic** and must be handled with utmost care.

**We recommend that:**

- Inhalation is avoided.
- Protection is used to avoid contact with skin and eyes.

**To avoid exposure of persons handling shipments:**

- Check that containers are adequately closed.
- Pack containers in pressurebag or another suitable bag.

**Technical information**

Formaldehyde 4% aqueous solution, buffered Q path® Formalin 10%

**Hazard pictograms**



Signal word: **Danger**

**Hazard statements**

H350 May cause cancer

H341 Suspected of causing genetic defects

H302 Harmful if swallowed

H312 Harmful in contact with skin

H332 Harmful if inhaled

H317 May cause an allergic skin reaction